

Silver Nitrate, Volumetric Solution, 0.171N

C7620	1L	67.08
C7620	4L	177.17
C7620	20L	730.40

Traceable to NIST standards

Silver Nitrate, Volumetric Solution, 0.1 N

79692-320	500 mL	34.34
79692-360	1 L	60.95
79692-540	4 L	158.47
79692-700	20 L	451.57

1mL = 16.99mg AgNO₃

Normality0.0995 - 0.1005

Traceable to NIST standards

Silver Nitrate, Volumetric Solution, 0.05 N

79695-320	500 mL	24.77
79695-360	1 L	34.59
79695-540	4 L	96.07

1mL = 8.494mg AgNO₃

Normality0.0498 - 0.0502

Traceable to NIST standards

Silver Nitrate, 0.0282 N

79697-320	500 mL	19.29
79697-360	1 L	25.01
79697-540	4 L	60.07

1 mL = 1 mg Cl

Normality0.0281 - 0.0283

Traceable to NIST standards

Silver Nitrate, Volumetric Solution, 0.0192 N

For cyanide

79703-320	500 mL	17.53
79703-360	1 L	26.07

1 mL = 2 mg CN

APHA 4500-CN⁻ D

Normality0.0187 - 0.0197

Traceable to NIST standards

Silver Nitrate, Volumetric Solution, 0.0171 N

79676-360	1 L	26.07
79676-540	4 L	62.52
79676-700	20 L	182.78

1 mL = 1 mg NaCl

Normality0.0166 - 0.0176

Traceable to NIST standards

Silver Nitrate, Volumetric Solution, 0.0141 N

79705-320	500 mL	14.96
79705-360	1 L	19.29
79705-540	4 L	46.40

1 mL = 0.5 mg Cl

Normality0.0136 - 0.0146

Traceable to NIST standards

Silver Nitrate, 10% Solution, W/V

79662-320	500 mL	137.56
-----------	--------	--------

Silver Nitrate, 5% Solution, W/V

79665-320	500 mL	112.19
-----------	--------	--------