

Fehling Solution A

Copper solution

40640-360	1 L	20.80
40640-540	4 L	45.11
40640-700	20 L	134.51

Fehling Solution B

Alkaline solution

40643-360	1 L	29.21
40643-540	4 L	73.86
40643-700	20 L	190.49

Ferric Ammonium Sulfate, Reagent, A.C.S.

$\text{NH}_4\text{Fe}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$ FW 482.20
CAS 7783-83-7

40710-300	500 g	83.42	6X500 g	400.47
-----------	-------	-------	---------	--------

A.C.S. Specifications

Appearance	Pale violet crystals
Assay $[\text{NH}_4\text{Fe}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}]$	98.5 - 102.0%
Insoluble matter	0.01% max.
Chloride (Cl)	0.001% max.
Nitrate (NO_3)	0.01% max.
Copper (Cu)	0.003% max.
Zinc (Zn)	0.003% max.
Ferrous iron (Fe^{++})	P.T.
Substances not precipitated by NH_4OH	0.05% max.

Ferric Ammonium Sulfate, 5% Solution, W/V

40733-360	1 L	49.09
-----------	-----	-------

Ferric Ammonium Sulfate, 6% Solution, W/V

40730-360	1 L	49.09
40730-540	4 L	72.46

Ferric Chloride, Anhydrous, Reagent

FeCl_3				FW 162.22
CAS 7705-08-0				
40790-140	100 g	28.17	6X100 g	135.34
40790-300	500 g	87.42	6X500 g	419.63
40790-580	10 kg	257.10		

Specifications

Assay (FeCl_3)	98.0% min.
Appearance	P.T.
Insoluble in HCl	0.5% max.
Ferrous Chloride	0.5% max.
Manganese (Mn)	0.5% max.

Ferric Chloride, 0.025% Solution

A.P.H.A. Method 5210 B.(BOD)

R2335	500 ml	19.05
R2335	1 L	28.16

Ferric Chloride 40% W/V Solution

C2143	1L	80.98
C2143	4L	254.64